

Human GLB1/ELNR1 (His tag) recombinant protein



Catalog Number: 518103

General Information

Gene Name Synonym

Beta-Galactosidase, Acid Beta-Galactosidase, Lactase, Elastin Receptor 1, GLB1, ELNR1

Protein Construction

A DNA sequence encoding the Human GLB1/ELNR1 Leu24-Val677 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human GLB1/ELNR1 (His tag) recombinant protein