

# H3N2 HA (His Tag) recombinant protein



Catalog Number: 501503

## General Information

### Protein Construction

A DNA sequence encoding the Influenza A virus (A/Hong Kong/CUHK31987/2011(H3N2)) hemagglutinin (AGC13545.1) (Met1-Trp530), termed as HA, was expressed with a polyhistidine tag at the C-terminus.

### Organism

H3N2

### Expression Host

Baculovirus-Insect Cells

## QC Testing

### Activity

Measured by its ability to agglutinate human red blood cells.

HA titer is 0.8-4  $\mu\text{g}/\text{mL}$  for 1%GRBC

### Purity

> 90 % as determined by SDS-PAGE.

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Gln 17

### Molecular Mass

The recombinant hemagglutinin of Influenza A virus (A/Hong Kong/CUHK31987/2011(H3N2)) consists 525 amino acids and predicts a molecular mass of 59.3 kDa.

### Formulation

Lyophilized from sterile 20 mM Tris, 500 mM NaCl, 10 % glycerol, pH 7.0.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Influenza A H3N2 (A/Hong Kong/CUHK31987/2011) Hemagglutinin / HA Protein (His Tag)