

HNRNPC antibody



Catalog Number: 111503

Product name

HNRNPC antibody

WB: 1:1000-1:10000

IP: 1:200-1:2000

Specificity

Human, Mouse, Rat; other species not tested.


IHC: 1:20-1:200

IF: 1:10-1:100

Antibody description

HNRNPC Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, human brain tissue, K-562 cells. Positive FC detected in HeLa cells. Positive IHC detected in human lymphoma tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 32 kDa

Validations




Preparation

This antibody was obtained by immunization of HNRNPC recombinant protein (Accession Number: XM_011536712). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111503(HNRNPC antibody) at dilution of 1:500

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:111503(HNRNPC antibody) at dilution of 1:100 (under 10x lens)

Applications


ELISA, WB, IHC, IF, FC, IP

Dilutions

Recommended Dilution:


HNRNPC antibody

Catalog Number: 111503




Immunofluorescent analysis of HeLa cells, using HNRNPC antibody Catalog No:111503 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit

IP Result of anti-HNRNPC (IP:Catalog No:111503, 4ug; Detection:Catalog No:111503 1:500) with HeLa cells lysate 1080ug.



1X10⁶ HeLa cells were stained with 0.2ug HNRNPC antibody (Catalog No:111503, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.



IgG (red).