

Rat CD68 (Fc Tag) recombinant protein



Catalog Number: 502503

General Information

Gene Name Synonym

Protein Cd68

Protein Construction

A DNA sequence encoding the rat CD68 (Q4FZY1) (Met1-Ser295) was expressed, fused with the Fc region of human IgG1 at the C-terminus.

Organism

Rat

Expression Host

Human Cells

QC Testing

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μ g of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 21

Molecular Mass

The recombinant rat CD68/Fc is a disulfide-linked

homodimer. The reduced monomer comprises 517 amino acids and has a predicted molecular mass of 56.5 kDa. The apparent molecular mass of the protein is approximately 87-117 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Rat CD68 / Macrosialin Protein (Fc Tag) SDS-PAGE