

Rac1 antibody



Catalog Number: 114503

Product name

Rac1 antibody

Specificity

Human; other species not tested.

Antibody description

Rac1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 20-21 kDa

Preparation

This antibody was obtained by immunization of Rac1 recombinant protein (Accession Number: NM_006908). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114503(RAC1 antibody) at dilution of 1:1000