

GCHFR antibody

Catalog Number: 110903

Product name

GCHFR antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

GCHFR Rabbit Polyclonal antibody. Positive IHC detected in human cervix. Positive WB detected in mouse skin tissue. Observed molecular weight by Western-blot: 10kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_005258). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



mouse skin tissue were subjected to SDS PAGE followed by western blot with Catalog No:110903(GCHFR antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human cervix using Catalog No:110903(GCHFR antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human cervix using Catalog No:110903(GCHFR antibody) at dilution of 1:100 (under 40x lens)