

Mouse CD99L2 (His Tag) recombinant protein



Catalog Number: 503903

General Information

Gene Name Synonym

MIC2-like protein 1

Protein Construction

A DNA sequence encoding the mouse CD99L2 (NP_612182.1) extracellular domain (Met 1-Ala 164) was expressed, fused with a polyhistidine tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Activity

Measured by its ability to bind biotinylated recombinant mouse CD99L2 in functional ELISA.

Purity

> 98 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 26

Molecular Mass

The secreted recombinant mouse CD99L2 comprises 150 amino acids and has a predicted molecular mass of 16.3 kDa. As a result of glycosylation, the apparent molecular mass of rmCD99L2 is approximately 30-40 kDa in SDS-PAGE under reducing conditions.

Formulation

- Lyophilized from sterile PBS, pH 7.4
1. 5 % trehalose and mannitol are added as protectants before lyophilization.
 2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse CD99L2 / MIC2L1 Protein (His Tag) SDS-PAGE