

PIK3C2G antibody



Catalog Number: 113903

Product name

PIK3C2G antibody

Immunogen

[Human PIK3C2G Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

PIK3C2G Rabbit Polyclonal antibody. Positive WB detected in A549 cells. Observed molecular weight by Western-blot: 150 kDa

Preparation

This antibody was obtained by immunization of PIK3C2G recombinant protein (Accession Number: XM_047429009). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:113903(PIK3C2G Antibody) at dilution of 1:600