

RPL9 antibody

Catalog Number: 114903

Product name

RPL9 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

RPL9 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 22kd

Preparation

This antibody was obtained by immunization of RPL9 recombinant protein (Accession Number: NM_000661). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114903(RPL9 antibody) at dilution of 1:800