

DGAT1 antibody

Catalog Number: 109903

Product name

DGAT1 antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DGAT1 Rabbit Polyclonal antibody. Positive IP detected in mouse heart tissue. Positive WB detected in mouse heart tissue. Positive IHC detected in human hepatocirrhosis tissue, human pancreas cancer tissue. Observed molecular weight by Western-blot: 50-57 kDa

Preparation

This antibody was obtained by immunization of DGAT1 recombinant protein (Accession Number: NM_012079). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:109903(DGAT1 Antibody) at dilution of 1:1000



IP Result of anti-DGAT1 (IP:Catalog No:109903, 3ug; Detection:Catalog No:109903 1:300) with mouse heart tissue lysate 4000ug.




Immunohistochemistry of paraffin-embedded human hepatocirrhosis tissue slide using Catalog

DGAT1 antibody

Catalog Number: 109903

No:109903(DGAT1 Antibody) at dilution of 1:50
(under 10x lens)



Immunohistochemistry of paraffin-embedded
human hepatocirrhosis tissue slide using Catalog
No:109903(DGAT1 Antibody) at dilution of 1:50
(under 40x lens)