

Anti-OR52E6 antibody

Catalog Number: 176203

Product name

Anti-OR52E6 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to OR52E6

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or52e6. aa range:264-313

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations



Fig1.; Western Blot analysis of various cells using Olfactory receptor 52E6 Polyclonal Antibody diluted at 1: 1000