

# Anti-Histone H2B (7E2) Mouse antibody



Catalog Number: 168203

## Product name

Anti-Histone H2B (7E2) Mouse antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Histone H2B (7E2) Mouse monoclonal antibody

## Preparation

Antigen: Synthetic Peptide of Histone H2B

## Formulation

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


WB, IHC-P

## Dilutions


WB: 1/1000-3000

IF: 1/200

## Validations



Western blot detection of Histone H2B in HeLa and COS7 cell lysates using Histone H2B mouse mAb (1:100000 diluted). Predicted band size: 14kDa. Observed band size: 14kDa.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma using Histone H2B mouse mAb diluted at 1:500.