

# Mouse TNFR2/CD120b/TNFRSF1B (Fc Tag) recombinant protein



Catalog Number: 501603

## General Information

### Gene Name Synonym

Tumor necrosis factor receptor 2; Tumor necrosis factor receptor type II; p75; p80 TNF-alpha receptor

### Protein Construction

A DNA sequence encoding the extracellular domain of mouse TNFRSF1B (NP\_035740.2) (Met 1-Gly 258) was fused with the Fc region of human IgG1 at the C-terminus.

### Organism

Mouse

### Expression Host

Human Cells

## QC Testing

### Activity

Measured in a cytotoxicity assay using L929 mouse fibrosarcoma cells in the presence of the metabolic inhibitor actinomycin D. The ED<sub>50</sub> for this effect is typically 4-16 ng/mL in the presence of 0.5 ng/mL of Recombinant Mouse TNF $\alpha$ .

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu$ g of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Val 23

### Molecular Mass

The recombinant mouse TNFRSF1B/Fc is a disulfide-linked homodimer. The reduced monomer consists of 477 amino acids and has a predicted molecular mass of 52.3 kDa. In SDS-PAGE under reducing conditions, rm TNFRSF1B/Fc monomer migrates as an approximately 65 kDa band due to glycosylation.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Mouse TNFR2 / CD120b / TNFRSF1B Protein (Fc

**Mouse TNFR2/CD120b/TNFRSF1B (Fc Tag)  
recombinant protein**



Catalog Number: 501603

---

Tag) SDS-PAGE