

BAF250 Polyclonal Antibody



Catalog Number: 163303

Product name

BAF250 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to BAF250

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 1300-1400 of mouse BAF250 (NP_001074288.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - BAF250 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using BAF250 antibody at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 15s.