

CHTF18 antibody

Catalog Number: 109303

Product name

CHTF18 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

CHTF18 Rabbit Polyclonal antibody. Positive IHC detected in human medulloblastoma. Positive WB detected in Y79 cells. Observed molecular weight by Western-blot: 110kd

Preparation

This antibody was obtained by immunization of CHTF18 recombinant protein (Accession Number: NM_022092). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:109303(CHTF18 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human medulloblastoma using Catalog No:109303(CHTF18 antibody) at dilution of 1:100 (under 10x lens)