

HSPBP1 antibody

Catalog Number: 111703

Product name

HSPBP1 antibody

ELISA, WB, IHC

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IHC: 1:20-1:200


Specificity

human, mouse; other species not tested.

Antibody description

HSPBP1 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer. Positive WB detected in mouse testis tissue, human brain tissue, human heart tissue, human testis tissue. Observed molecular weight by Western-blot: 40kd

Validations




Preparation

This antibody was obtained by immunization of HSPBP1 recombinant protein (Accession Number: NM_012267). Purification method: Antigen affinity purified.

mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:111703(HSPBP1 antibody) at dilution of 1:500

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:111703(HSPBP1 antibody) at dilution of 1:100 (under 25x lens)

Applications