

PCF11 antibody



Catalog Number: 113703

Product name

PCF11 antibody

Immunogen

[Human PCF11 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

PCF11 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, HeLa cells. Observed molecular weight by Western-blot: 173 kDa

Preparation

This antibody was obtained by immunization of PCF11 recombinant protein (Accession Number: NM_015885). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:113703(PCF11 antibody) at dilution of 1:1000