

Anti-H2AFM/HIST1H2AB antibody



Catalog Number: 104703

Product name

Anti-H2AFM/HIST1H2AB antibody

Specificity

Human H2AFM / HIST1H2AB

Antibody description

Rabbit polyclonal to H2AFM/HIST1H2AB

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human H2AFM / HIST1H2AB, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1-10 µg/ml

Validations



Lanes	A	B	C
Sample (whole cell lysate)	K562	MDA-MB-231	HEK293
Sample Volume (µg/lane)	30	30	30
Gel	13% SDS-PAGE reducing gel		
Recommended Concentration	1-10 µg/ml		
Secondary Antibody	Dylight 800 Labeled Antibody to Rabbit IgG (H+), at 1:5000 dilution.		
Explanation	Developed using Odyssey imaging system. Predicted band size: 14 kDa Observed band size: 14 kDa Additional bands: at ~22 kDa (We are unsure as to the identity of these extra bands.)		

H2AFM / HIST1H2AB Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot