

## General Information

### Protein Construction

A DNA sequence encoding the rat IFNA5 (XP\_001076062.1) (Met1-Glu189) was expressed, fused with a polyhistidine tag at the C-terminus.

### Organism

Rat

### Expression Host

Human Cells

## QC Testing

### Activity

Measured in antiviral assays using L929 cells infected with vesicular stomatitisvirus (VSV). The ED<sub>50</sub> for this effect is typically 10-50 pg/mL.

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Cys 24

### Molecular Mass

The recombinant rat IFNA5 comprises 179 amino acids and predicts a molecular mass of 20.9 kDa. The apparent molecular mass of the recombinant protein is approximately 21 kDa in SDS-PAGE under reducing conditions due to glycosylation.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Rat IFNA5 / IFNaG Protein (His Tag) SDS-PAGE