

# STEAP1 antibody



Catalog Number: 115703

---

## Product name

STEAP1 antibody

## Specificity

human; other species not tested.

## Antibody description

STEAP1 Rabbit Polyclonal antibody. Positive WB detected in A431 cells, DU 145 cells. Observed molecular weight by Western-blot: 38-40kd

## Preparation

This antibody was obtained by immunization of STEAP1 recombinant protein (Accession Number: NM\_012449). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:115703(STEAP1 antibody) at dilution of 1:500