

Anti-PPCEL/PREPL Rabbit antibody



Catalog Number: 167703

Product name

Anti-PPCEL/PREPL Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

PPCEL/PREPL Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human PPCEL/PREPL

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/1000

ICC/IF: 1/200

FC: 1/20

IP: 1/20

Validations



Western blot detection of PPCEL/PREPL in HeLa,C6,Ramos using PPCEL/PREPL Rabbit mAb(1:1000 diluted)