

Anti-TMEM132B antibody

Catalog Number: 175013

Product name

Polyclonal

Anti-TMEM132B antibody

Ig Type**Specificity**

Rabbit IgG

Human, Mouse

Applications**Antibody description**

WB

Rabbit polyclonal antibody to TMEM132B


Dilutions**Preparation**

WB:1:500□1:1,000

This antigen of this antibody was synthetic peptide within n-terminal human tmem132b.

Validations**Formulation**

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

**Storage**

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

Fig1: Western blot analysis on D3 using anti-TMEM132B polyclonal antibody.

Clonality