

Anti-Phospho-CBL (Tyr774) Rabbit antibody



Catalog Number: 169013

Product name

Monoclonal

Anti-Phospho-CBL (Tyr774) Rabbit antibody

Ig Type

Specificity

Rabbit IgG

Human

Applications

Antibody description

WB, FC, IP

Phospho-CBL (Tyr774) Rabbit monoclonal antibody

Dilutions

WB: 1/1000

Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Tyr774 of human CBL

Formulation


Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Validations



Western blot detection of Phospho-CBL (Tyr774) in K562,C6,HeLa cell lysates using Phospho-CBL (Tyr774) Rabbit mAb(1:1000 diluted).Predicted band size:100kDa.Observed band size:120kDa.