

# CD24 antibody



Catalog Number: 109013

---

## Product name

CD24 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CD24 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, HeLa cells, rat spleen tissue. Observed molecular weight by Western-blot: 30-70 kDa

## Preparation

This antibody was obtained by immunization of CD24 recombinant protein (Accession Number: BC007674). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109013(CD24 antibody) at dilution of 1:1000