

Anti-Cytokeratin 7 (6E10) Mouse antibody



Catalog Number: 168413

Product name

Anti-Cytokeratin 7 (6E10) Mouse antibody

WB: 1/500 - 1/2000

FC: 1/200 - 1/400

Specificity

Human

ELISA: 1/10000

Antibody description

Cytokeratin 7 (6E10) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human CK7 expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type


Mouse IgG1

Applications


WB, FC

Dilutions

Validations



Flow cytometric analysis of HeLa cells using anti-CK7 mAb (green) and negative control (purple).



Western blot detection of Keratin 7 in A549 and HeLa cell lysates using Keratin 7 mouse mAb (1:500 diluted). Predicted band size: 51 kDa. Observed band size: 51 kDa.