

Human STAT5B (His tag) recombinant protein



Catalog Number: 519413

General Information

Gene Name Synonym

Signal Transducer and Activator of Transcription 5B, STAT5B

Protein Construction

A DNA sequence encoding the Human STAT5B Met1-Thr321 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of PBS, 50% Glycerol, 1mM DTT, pH 7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human STAT5B (His tag) recombinant protein