

Elsilimomab (Synonyms: B-E8; MAB B-E8)



Catalog Number: 176813

Product name

Elsilimomab (Synonyms: B-E8; MAB B-E8)

Specificity

Human IL-6 / IFN β 2

Antibody description

Anti-IL-6 / IFN β 2 Antibody (Elsilimomab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG1

Applications


ELISA

Validations




SDS-PAGE

Anti-IL-6 / IFN β 2 Antibody (elsilimomab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-IL-6 / IFN β 2 Antibody (elsilimomab) is 98.96%, determined by SEC-



HPLC.

Bioactivity: ELISA

Immobilized human IL-6 His at 2 μ g/mL can bind Anti-IL-6 / IFN β 2 Antibody (elsilimomab), EC $_{50}$ =0.0162 μ g/mL.