

# Anti-NF1B antibody



Catalog Number: 101113

## Product name

Anti-NF1B antibody

## Specificity

human NF1B

## Antibody description

Rabbit polyclonal to NF1B

## Preparation

Produced in rabbits immunized with A synthetic peptide corresponding to the C-terminus of the human NF1B, and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 5-10 µg/ml

## Validations



Items	Lanes	A
Sample (whole cell lysate)		293T
Sample Volume (µg/lane)		30
Gel		13% SDS PAGE reducing gel
Recommended Concentration		5-10 µg/ml
Secondary Antibody		Dylight 800-labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution.
Developed using Odyssey imaging system.		
Explanation		Predicted band size : 47 kDa Observed band size : 47 kDa

NF1B Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot