

ALOX5 antibody



Catalog Number: 107113

Product name

ALOX5 antibody

Specificity

Human, Pig; other species not tested.

Antibody description

ALOX5 Mouse Monoclonal antibody. Positive WB detected in pig spleen tissue. Observed molecular weight by Western-blot: 78 kDa

Preparation

This antibody was obtained by immunization of ALOX5 recombinant protein (Accession Number: NM_000698). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



pig spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:107573(ALOX5 Antibody) at dilution of 1:2000