

# CDCA2 antibody

Catalog Number: 109113

## Product name

CDCA2 antibody

## Specificity

Human; other species not tested.

## Antibody description

CDCA2 Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells. Positive WB detected in Jurkat cells, HepG2 cells, MCF7 cells. Observed molecular weight by Western-blot: 125 kDa

## Preparation

This antibody was obtained by immunization of CDCA2 recombinant protein (Accession Number: BC036214). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109113(CDCA2 antibody) at dilution of 1:500



IP Result of anti-CDCA2 (IP:Catalog No:109113, 4ug; Detection:Catalog No:109113 1:800) with Jurkat cells lysate 2800ug.