

# OAF Polyclonal Antibody



Catalog Number: 161513

## Product name

OAF Polyclonal Antibody

## Specificity

Mouse, Rat

## Antibody description

Polyclonal antibody to OAF

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 28-150 of human OAF (NP\_848602.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



## Western blot - OAF Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using OAF antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.