

HNRNPU antibody



Catalog Number: 111513

Product name

HNRNPU antibody

WB: 1:1000-1:10000

IP: 1:200-1:2000

Specificity


Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

Antibody description

HNRNPU Rabbit Polyclonal antibody. Positive FC detected in HepG2 cells. Positive IHC detected in human liver tissue. Positive WB detected in HEK-293 cells, HeLa cells, HepG2 cells, Jurkat cells, K-562 cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 120 kDa

Validations




Preparation

This antibody was obtained by immunization of HNRNPU recombinant protein (Accession Number: XM_017001117). Purification method: Antigen affinity purified.

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111513(HNRNPU antibody) at dilution of

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

1:500

Clonality

Polyclonal

IP Result of anti-HNRNPU (IP:Catalog No:111513, 3ug; Detection:Catalog No:111513 1:800) with HeLa cells lysate 2000ug.

Ig Type


Rabbit IgG

Applications

ELISA, WB, IP, IHC, FC

Dilutions

Recommended Dilution:




HNRNPU antibody



Catalog Number: 111513

Immunohistochemical of paraffin-embedded human liver using Catalog No:111513(HNRNPU antibody) at dilution of 1:100 (under 40x lens)



1X10⁶ HepG2 cells were stained with 0.2ug HNRNPU antibody (Catalog No:111513, red) and control antibody (blue). Fixed with 4% PFA blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.