

# UBE1 antibody

Catalog Number: 116513

## Product name

UBE1 antibody

IHC: 1:20-1:200

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

UBE1 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive WB detected in HL-60 cells, HEK-293 cells, HeLa cells, mouse colon tissue, mouse spleen tissue, rat spleen tissue. Observed molecular weight by Western-blot: 118 kDa

## Preparation

This antibody was obtained by immunization of UBE1 recombinant protein (Accession Number: XM\_017029781). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


## Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:116513(UBE1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:116513(UBE1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human colon cancer using Catalog

# UBE1 antibody



Catalog Number: 116513

---

No:116513(UBE1 antibody) at dilution of 1:50 (under 40x lens)