

RSPO3 antibody



Catalog Number: 107513

Product name

RSPO3 antibody

Specificity

Human; other species not tested.

Antibody description

RSPO3 Mouse Monoclonal antibody. Positive WB detected in human skeletal muscle tissue, A549 cells. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of RSPO3 recombinant protein (Accession Number: NM_032784). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:107513(RSPO3 Antibody) at dilution of 1:2000