

Anti-ELK1 Mouse antibody



Catalog Number: 168513

Product name

Anti-ELK1 Mouse antibody

Specificity

Human

Antibody description

ELK1 Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of ELK1 expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


WB

Dilutions

WB: 1/500 - 1/2000

ELISA: 1/10000

Validations



Western blot analysis using ELK1 mouse mAb against truncated ELK1 recombinant protein (1) and K562 cell lysate (2).