

TIMM10B Polyclonal Antibody



Catalog Number: 165913

Product name

TIMM10B Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to TIMM10B

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-103 of human TIMM10B (NP_036324.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100

Validations



Western blot - TIMM10B Polyclonal Antibody