

# TCN2 antibody

Catalog Number: 115913

## Product name

TCN2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TCN2 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue. Positive WB detected in A375 cells, mouse small intestine tissue, rat liver tissue, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 48kd

## Preparation

This antibody was obtained by immunization of TCN2 recombinant protein (Accession Number: NM\_001184726). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:115913(TCN2 antibody) at dilution of 1:100



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:115913(TCN2 antibody) at dilution of 1:50 (under 10x lens)