

# DGKK-Specific antibody



Catalog Number: 109913

## Product name

DGKK-Specific antibody

## Specificity

human; other species not tested.

## Antibody description

DGKK-Specific Rabbit Polyclonal antibody. Positive IF detected in HepG2. Positive WB detected in human testis tissue. Observed molecular weight by Western-blot: 135-140kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_001013742). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109913(DGKK-Specific antibody) at dilution of



Immunofluorescent analysis of HepG2 cells, using DGKK antibody Catalog No:109913 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).