

TPPP2 antibody

Catalog Number: 116213

Product name

TPPP2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

TPPP2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 18kd

Preparation

This antibody was obtained by immunization of TPPP2 recombinant protein (Accession Number: XM_017020966). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116213(TPPP2 antibody) at dilution of 1:100