

Anti-CD79a Rabbit antibody



Catalog Number: 167213

Product name

Anti-CD79a Rabbit antibody

Specificity

Human

Antibody description

CD79a Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human CD79a

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/1000-1/5000

IHC: 1/50-1/200

ICC/IF: 1/50-1/200

FC: 1/30

Validations



Western blot analysis of CD79a expression in human spleen lysate.