

ASB5 antibody



Catalog Number: 108213

Product name

ASB5 antibody

Immunogen

[Human ASB5 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

ASB5 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, mouse heart tissue. Observed molecular weight by Western-blot: 36-40 kDa

Preparation

This antibody was obtained by immunization of ASB5 recombinant protein (Accession Number: XM_005262759). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:108213(ASB5 Antibody) at dilution of 1:300