

# Human CXCL8/IL8 recombinant protein



Catalog Number: 518213

## General Information

### Gene Name Synonym

Interleukin-8, IL-8, C-X-C Motif Chemokine 8, Emotakin, Granulocyte Chemotactic Protein 1, GCP-1, Monocyte-Derived Neutrophil Chemotactic Factor, MDNCF, Monocyte-Derived Neutrophil-Activating Peptide, MONAP, Neutrophil-Activating Protein 1, NAP-1, Protein 3-10C, T-Cell Chemotactic Factor, IL8, CXCL8

### Protein Construction

A DNA sequence encoding the Human CXCL8/IL8 Ala23-Ser99 is expressed.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Store for up to 12 months at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  as lyophilized powder. Aliquots of reconstituted

samples are stable at  $< -20^{\circ}\text{C}$  for 3 months.

### Formulation

Lyophilized from a  $0.2 \mu\text{m}$  filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

## Usage Guide


### Storage

Lyophilized protein should be stored at  $< -20^{\circ}\text{C}$ , though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at  $4-7^{\circ}\text{C}$  for 2-7 days.

### Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than  $100 \mu\text{g/ml}$ . Dissolve the lyophilized protein in ddH<sub>2</sub>O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

## SDS-PAGE



Human CXCL8/IL8 recombinant protein