

General Information

Gene Name Synonym

N-Sulphoglucosamine Sulphohydrolase,
Sulfoglucosamine Sulfamidase, Sulphamidase,
SGSH, HSS

Protein Construction

A DNA sequence encoding the Human SGSH/HSS Arg21-Leu502 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 1mM GaCl₂, 10% Glycerol, pH 7.5.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human SGSH/HSS (His tag) recombinant protein