

Mouse TIMP-1/TIMP1 recombinant protein



Catalog Number: 501613

General Information

Gene Name Synonym

Collagenase inhibitor 16C8 fibroblast; Erythroid-potentiating activity; TPA-S1; TPA-induced protein; Tissue inhibitor of metalloproteinases 1

Protein Construction

A DNA sequence encoding the mouse TIMP1 (NP_035723.2) (Met 1-Arg 205) was expressed, with a C-terminal polyhistidine tag.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Activity

Measured by its ability to inhibit human MMP-2 cleavage of a fluorogenic peptide substrate MCA-PLGL-DPA-AR-NH₂(R&D Systems, Catalog # ES001). The IC₅₀ value is < 6 nM.

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Cys 25

Molecular Mass

The secreted recombinant mouse TIMP1 consists of 181 amino acids and has a calculated molecular mass of 21 kDa. It migrates as an approximately 30 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse TIMP1 / TIMP / CLGI Protein SDS-PAGE