

MID1IP1 antibody



Catalog Number: 112613

Product name

MID1IP1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MID1IP1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 23 kDa, 46 kDa

Preparation

This antibody was obtained by immunization of MID1IP1 recombinant protein (Accession Number: NM_001098790). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112613(MID1IP1 Antibody) at dilution of 1:300