

PRKCG Polyclonal Antibody



Catalog Number: 165613

Product name

PRKCG Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to PRKCG

Preparation

Antigen: A synthetic peptide of human PRKCG

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB, IF

Dilutions


WB 1:500 - 1:2000
IF 1:50 - 1:200

Validations



Western blot - PRKCG Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using PRKCG antibody at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time:



Immunofluorescence - PRKCG Polyclonal Antibody

Immunofluorescence analysis of rat brain cells using PRKCG antibody at dilution of 1:100. Blue: DAPI for nuclear staining.