

# KCNA1 Polyclonal Antibody



Catalog Number: 166613

---

## Product name

KCNA1 Polyclonal Antibody

## Specificity

Mouse

## Antibody description

Polyclonal antibody to KCNA1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 426-495 of human KCNA1 (NP\_000208.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IF

## Dilutions

IF 1:50 - 1:200

## Validations



Immunofluorescence - KCNA1 Polyclonal Antibody

Immunofluorescence analysis of mouse brain using KCNA1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.