

Anti-CCL2 antibody



Catalog Number: 176613

Product name

Rabbit IgG

Anti-CCL2 antibody

Applications

Specificity

WB

Mouse, Human

Dilutions

Antibody description

WB: 1:500-1:1,000

Rabbit polyclonal antibody to CCL2

Validations

Preparation

This antigen of this antibody was synthetic peptide within human ccl2 aa 30-70 / 99.

Formulation

Liquid, 1*TBS (pH7.4), 0.2% BSA, 50% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type




Fig1:: Western blot analysis of CCL2 on mouse testis tissue lysate. Proteins were transferred to a PVDF membrane and blocked with 5% NFDM/TBST for 1 hour at room temperature. The primary antibody (1/500) was used in 5% NFDM/TBST at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody (HA1001) at 1:200,000 dilution was used for 1 hour at room temperature.