

# CST6 antibody

Catalog Number: 109613

---

## Product name

CST6 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CST6 Rabbit Polyclonal antibody. Positive IHC detected in mouse liver.

## Preparation

This antibody was obtained by immunization of CST6 recombinant protein (Accession Number: NM\_001323). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded mouse liver using Catalog No:109613(CST6 antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemical of paraffin-embedded mouse liver using Catalog No:109613(CST6 antibody) at dilution of 1:50 (under 10x lens)