

# Anti-SLC17A2 antibody



Catalog Number: 176313

---

**Product name**

Polyclonal

Anti-SLC17A2 antibody

**Ig Type****Specificity**

Rabbit IgG

Human, Mouse

**Applications****Antibody description**

WB

Rabbit polyclonal antibody to SLC17A2

**Dilutions****Preparation**

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human slc17a2. aa range:213-262

**Validations****Formulation**

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

**Storage**

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.


**Clonality**

Fig1:: Western Blot analysis of various cells using SLC17A2 Polyclonal Antibody